

SURF1 antibody

Catalog Number: 115746

Product name

SURF1 antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SURF1 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Positive IHC detected in human heart tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 30-31 kDa

Preparation

This antibody was obtained by immunization of SURF1 recombinant protein (Accession Number: NM_003172). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115746(SURF1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human heart using Catalog No:115746(SURF1 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human heart using Catalog No:115746(SURF1

SURF1 antibody

Catalog Number: 115746

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:115746(SURF1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)