

# Anti-HPRT Rabbit antibody



Catalog Number: 166746

## Product name

Anti-HPRT Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

HPRT Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human HPRT

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions


WB: 1/2000-1/10000

IP: 1/20

## Validations



Western blot detection of HPRT in HeLa, CHO-K1, C6, Rat Brain, Jurkat cell lysates using HPRT Rabbit pAb(1:500 diluted). Predicted band size: 25kDa. Observed band size: 25kDa.



Western blot detection of HPRT in K562, 3T3, HeLa cell lysates using HPRT Rabbit pAb(1:1000 diluted). Predicted band size: 25kDa. Observed band size: 25kDa.