

NDOR1 antibody

Catalog Number: 113056

Product name

NDOR1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NDOR1 Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue, human liver tissue, human lung tissue, human skin tissue, mouse testis tissue, mouse thymus tissue. Positive IHC detected in human lymphoma. Observed molecular weight by Western-blot: 62kd

Preparation

This antibody was obtained by immunization of NDOR1 recombinant protein (Accession Number: NM_014434). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:113056(NDOR1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:113056(NDOR1 antibody) at dilution of 1:50 (under 40x lens)