

Anti-HNF-4-alpha Rabbit antibody



Catalog Number: 168056

Product name

Rabbit IgG

Anti-HNF-4-alpha Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC, IP

Human

Dilutions

Antibody description

WB: 1/1000

HNF-4-alpha Rabbit monoclonal antibody

ICC/IF: 1/20

Preparation

FC: 1/100


Antigen: Recombinant protein of human HNF-4-alpha

IP: 1/20

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Validations



Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Western blot detection of HNF-4-alpha in A549, MCF-7 cell lysates using HNF-4-alpha Rabbit mAb(1:1000 diluted). Predicted band size: 53kDa. Observed band size: 53kDa.