

Anti-DCTN1/p150-glued (3D5) Mouse antibody



Catalog Number: 169056

Product name

WB, IP

Anti-DCTN1/p150-glued (3D5) Mouse antibody

Dilutions

Specificity

WB: 1/500

Human

Antibody description

DCTN1/p150-glued (3D5) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant Dynactin 1 protein fragments expressed in E.coli

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Mouse IgG2b

Applications

Validations



Immunoprecipitation analysis of HeLa cell lysates using Dynactin 1 mouse mAb.



Western blot detection of Dynactin 1 in K562, MCF7, 293T and HeLa cell lysates using Dynactin 1 mouse mAb (1:500 diluted). Predicted band size: 150KDa. Observed band size: 150KDa.