

# RAB6A antibody

Catalog Number: 114456

## Product name

RAB6A antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RAB6A Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 24 kDa

## Preparation

This antibody was obtained by immunization of RAB6A recombinant protein (Accession Number: NM\_002869). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114456(RAB6A antibody) at dilution of 1:500



IP Result of anti-RAB6A (IP:Catalog No:114456, 4ug; Detection:Catalog No:114456 1:500) with mouse brain tissue lysate 2640ug.