

# Phospho-ESPL1-S1126 pAb



Catalog Number: 166456

---

## Product name

Phospho-ESPL1-S1126 pAb

## Specificity

Human

## Antibody description

Polyclonal antibody to Phospho-ESPL1-S1126

## Preparation

Antigen: A synthetic phosphorylated peptide around S1126 of human ESPL1 (NP\_036423.4).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:20 - 1:100

## Validations



Immunofluorescence - Phospho-ESPL1-S1126 pAb

Immunofluorescence analysis of U2OS cells using Phospho-ESPL1-S1126 antibody . Blue: DAPI for nuclear staining.