

# Anti-MRPL35 antibody

Catalog Number: 176456

## Product name

Anti-MRPL35 antibody

## Specificity

Human, Mouse

## Antibody description

Rabbit polyclonal antibody to MRPL35

## Preparation

This antigen of this antibody was synthesized peptide derived from mrp-l35 . at aa range: 80-160

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1:; Western blot analysis of 293T COLO using MRP-L35 antibody. Antibody was diluted at 1:2000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000