

Anti-MRPL35 antibody



Catalog Number: 176456

Product name

Anti-MRPL35 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to MRPL35

Preparation

This antigen of this antibody was synthesized peptide derived from mrp-l35 . at aa range: 80-160

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1:; Western blot analysis of 293T COLO using MRP-L35 antibody. Antibody was diluted at 1:2000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000