

Anti-HIF-1 alpha Rabbit antibody



Catalog Number: 167456

Product name

Anti-HIF-1 alpha Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

HIF-1 alpha Rabbit polyclonal antibody

Preparation

Antigen: The antiserum was produced against synthesized peptide derived from human HIF-1alpha. AA range:328-377

Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ICC/IF, WB, IHC-P, ELISA

Dilutions


IF: 1/50-200

WB: 1/500 - 1/2000

IHC: 1/100 - 1/300

ELISA: 1/40000

Validations



Western Blot analysis of various cells using HIF-1u03b1 Polyclonal Antibody