

# P-cadherin antibody

Catalog Number: 107456

## Product name

P-cadherin antibody

## Specificity

Human; other species not tested.

## Antibody description

P-cadherin Mouse Monoclonal antibody. Positive WB detected in A431 cells. Positive IHC detected in human placenta tissue. Observed molecular weight by Western-blot: 120 kDa

## Preparation

This antibody was obtained by immunization of P-cadherin recombinant protein (Accession Number: NM\_001793). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:50-1:500


## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:107456(CDH3 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human placenta tissue slide using (CDH3 Antibody) at dilution of 1:500 (under 40x lens)



Immunohistochemistry of paraffin-embedded human placenta tissue slide using (CDH3 Antibody) at dilution of 1:500 (under 40x lens)