

# Anti-Cytochrome P450 17A1/CYP17A1 (6E1) Mouse antibody



Catalog Number: 168456

## Product name

Anti-Cytochrome P450 17A1/CYP17A1 (6E1)  
Mouse antibody

## Specificity

Human

## Antibody description

Cytochrome P450 17A1/CYP17A1 (6E1) Mouse  
monoclonal antibody

## Preparation

Antigen: Purified recombinant human Cytochrome  
P450 17A1 protein fragments expressed in E.coli.

## Formulation

Ascites containing with 0.03% Proclin300 and 50%  
glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of Cytochrome P450 17A1  
in HeLa, HepG2 and A549 cell lysates using  
Cytochrome P450 17A1 mouse mAb (1:1000  
diluted). Predicted band size: 60KDa. Observed  
band size: 60KDa.