

Mouse CD25/IL2RA (His Tag) recombinant protein



Catalog Number: 502856

General Information

Gene Name Synonym

p55

Protein Construction

A DNA sequence encoding the mouse IL2RA (NP_032393.3) extracellular domain (Met 1-Lys 236) was expressed, with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 98 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 22

Molecular Mass

The secreted recombinant mouse IL2RA comprises 226 amino acids with a predicted molecular mass of 26 kDa. As a result of glycosylation, the apparent molecular mass of the rm IL2RA is approximately 50-55 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse IL2RA / CD25 Protein (His Tag) SDS-PAGE