

# RTCD1 antibody

Catalog Number: 114856

---

## Product name

RTCD1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

RTCD1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 41-43kd

## Preparation

This antibody was obtained by immunization of RTCD1 recombinant protein (Accession Number: NM\_001130841). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114856(RTCD1 antibody) at dilution of 1:500