

Anti-PABPN1 Rabbit antibody



Catalog Number: 166856

Product name

Anti-PABPN1 Rabbit antibody

IHC: 1/20-1/100

ICC/IF: 1/20-1/50

Specificity

Human, Mouse

FC: 1/20

IP: 1/20

Antibody description

PABPN1 Rabbit polyclonal antibody

Validations

Preparation

Antigen: A synthetic peptide of human PABPN1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000-1/5000



Immunohistochemistry of PABPN1 in paraffin-embedded Human breast cancer tissue using PABPN1 Rabbit pAb at dilution 1/20



Western blot detection of PABPN1 in Jurkat, HeLa, COS7, Raji cell lysates using PABPN1 Rabbit pAb(1:1000 diluted). Predicted band size: 33KDa. Observed band size: 50KDa.