

Human CLIC3 (His tag) recombinant protein



Catalog Number: 518856

General Information

Gene Name Synonym

Chloride intracellular channel protein 3, CLIC3,

Protein Construction

A DNA sequence encoding the Human CLIC3 Met1-Arg236 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 10mM Tris, 0.1% Triton100, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human CLIC3 (His tag) recombinant protein