

Anti-COMMD9 antibody



Catalog Number: 100156

Product name

Anti-COMMD9 antibody

Immunogen

[Human COMMD9 \(His Tag\) recombinant protein](#)

Specificity

Human COMMD9

Antibody description

Rabbit monoclonal to COMMD9

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human COMMD9 (rh COMMD9; Q9P000-1; Met1-Lys198).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 0.1-3 µg/mL

Validations



COMMD9 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human COMMD9 in human lung with rabbit monoclonal antibody (0.2 µg/mL, formalin-fixed paraffin embedded



sections).

COMMD9 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human COMMD9 in human small intestine with rabbit monoclonal antibody (0.2 µg/mL, formalin-fixed paraffin embedded sections).