

GRAMD4 antibody



Catalog Number: 111156

Product name

Polyclonal

GRAMD4 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

GRAMD4 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 66, 140 kDa

Dilutions


Recommended Dilution:

Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of GRAMD4 recombinant protein (Accession Number: XM_047441209). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:111156(GRAMD4 Antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality