

ZNF287 antibody



Catalog Number: 117156

Product name

ZNF287 antibody

Specificity

Human; other species not tested.

Antibody description

ZNF287 Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 88-100 kDa

Preparation

This antibody was obtained by immunization of ZNF287 recombinant protein (Accession Number: NM_020653). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:117156(ZNF287 Antibody) at dilution of 1:300