

SPATC1 antibody



Catalog Number: 115556

Product name

SPATC1 antibody

Specificity

Human; other species not tested.

Antibody description

SPATC1 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of SPATC1 recombinant protein (Accession Number: BC053547). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115556(SPATC1 Antibody) at dilution of 1:600