

Mouse CD70/CD27L/TNFSF7 (Fc Tag) recombinant protein



Catalog Number: 500956

General Information

Gene Name Synonym

CD27 ligand; Tumor necrosis factor ligand superfamily member 7

Protein Construction

A DNA sequence encoding the mouse Cd70 (NP_035747.1) (Gln47-Pro159) was expressed with the Fc region of mouse IgG1 at the N-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> (12.5+86.3) % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp

Molecular Mass

The recombinant mouse Cd70 consists 385 amino acids and predicts a molecular mass of 43.1 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CD70 / CD27L / TNFSF7 Protein (Fc Tag)