

LGR6 Polyclonal Antibody



Catalog Number: 161956

Product name

LGR6 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to LGR6

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 700-800 of human LGR6 (NP_001017404.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - LGR6 Polyclonal Antibody

Western blot analysis of extracts of rat testis, using LGR6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.