

# ITPK1 antibody



Catalog Number: 111956

## Product name

ITPK1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ITPK1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HEK-293 cells, mouse eye tissue. Observed molecular weight by Western-blot: 45 kDa

## Preparation

This antibody was obtained by immunization of ITPK1 recombinant protein (Accession Number: NM\_001142594). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111956(ITPK1 antibody) at dilution of 1:300