

# PLEKHA3 antibody



Catalog Number: 113956

---

## Product name

PLEKHA3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PLEKHA3 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, Jurkat cells, K-562 cells, mouse brain tissue. Observed molecular weight by Western-blot: 34kd

## Preparation

This antibody was obtained by immunization of PLEKHA3 recombinant protein (Accession Number: NM\_019091). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:113956(PLEKHA3 antibody) at dilution of 1:200