

S100A12 antibody



Catalog Number: 114956

Product name

S100A12 antibody

Specificity

Human; other species not tested.

Antibody description

S100A12 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, human tonsillitis tissue.

Preparation

This antibody was obtained by immunization of S100A12 recombinant protein (Accession Number: NM_005621). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:114956(S100A12 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:114956(S100A12 Antibody) at dilution of 1:50 (under 40x lens)