

Anti-Proteasome subunit alpha type 6 Rabbit antibody



Catalog Number: 166956

Product name

Anti-Proteasome subunit alpha type 6 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Proteasome subunit alpha type 6 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Proteasome subunit alpha type 6

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/2000-1/10000

ICC/IF: 1/50

FC: 1/20-1/100

Validations



Western blot detection of Proteasome subunit alpha type 6 in K562,C6,3T3,HeLa cell lysates using Proteasome subunit alpha type 6 Rabbit pAb(1:1000 diluted).Predicted band size:27kDa.Observed band size:27kDa.