

# DNALI1 antibody



Catalog Number: 109956

## Product name

DNALI1 antibody

## Immunogen

[Human DNALI1 Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DNALI1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 30-33kd

## Preparation

This antibody was obtained by immunization of DNALI1 recombinant protein (Accession Number: NM\_003462). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:1000-1:10000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109956(DNALI1 antibody) at dilution of 1:3000