

HARS antibody

Catalog Number: 111256

Product name

HARS antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HARS Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HeLa cells, human kidney tissue, K-562 cells, mouse brain tissue, mouse heart tissue, mouse lung tissue. Positive IF detected in MCF-7. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of HARS recombinant protein (Accession Number: NM_002109). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111256(HARS antibody) at dilution of 1:500



Immunofluorescent analysis of MCF-7 cells, using HARS antibody Catalog No:111256 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).