

Product name

NMU antibody

Specificity

Human; other species not tested.

Antibody description

NMU Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human colon tissue, human small intestine tissue.

Preparation

This antibody was obtained by immunization of NMU recombinant protein (Accession Number: NM_006681). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:113256(NMU Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:113256(NMU Antibody) at dilution of 1:50 (under 40x lens)