

Torsin A antibody



Catalog Number: 116256

Product name

Torsin A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Torsin A Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in mouse brain tissue, HEK-293 cells, mouse kidney tissue, mouse liver tissue, rat brain tissue. Observed molecular weight by Western-blot: 42-50 kDa

Preparation

This antibody was obtained by immunization of Torsin A recombinant protein (Accession Number: NM_000113). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116256(TOR1A antibody) at dilution of 1:100



(under 25x lens)

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116256(TOR1A Antibody) at dilution of 1:1000