

# HBG1-Specific antibody



Catalog Number: 107256

## Product name

HBG1-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

HBG1-Specific Mouse Monoclonal antibody.  
Positive IHC detected in human placenta tissue.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000559).  
Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human placenta slide using (HBG1-Specific Antibody) at dilution of 1:200



Immunohistochemistry of paraffin-embedded human placenta slide using (HBG1-Specific Antibody) at dilution of 1:200