

Anti-PPP1A/PPP1CA (1H1) Mouse antibody



Catalog Number: 168256

Product name

Anti-PPP1A/PPP1CA (1H1) Mouse antibody

Specificity

Human

Antibody description

PPP1A/PPP1CA (1H1) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human PPP1A protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/50

Validations



Western blot detection of PPP1A in HepG2, MCF7 and HeLa cell lysates using PPP1A mouse mAb (1:1000 diluted). Predicted band size: 37KDa. Observed band size: 37KDa.



Immunocytochemistry of HeLa cells using anti-PPP1A mouse mAb diluted 1:50.