

# IL10RB antibody



Catalog Number: 111656

## Product name

IL10RB antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

IL10RB Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, Jurkat cells. Observed molecular weight by Western-blot: 42-45 kDa

## Preparation

This antibody was obtained by immunization of IL10RB recombinant protein (Accession Number: NM\_000628). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:111656(IL10RB antibody) at dilution of 1:300