

PDE12 antibody

Catalog Number: 113656

Product name

PDE12 antibody

IP: 1:200-1:2000

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

PDE12 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HEK-293 cells. Positive IP detected in mouse brain tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 65 kDa




Preparation

This antibody was obtained by immunization of PDE12 recombinant protein (Accession Number: BC034978). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113656(PDE12 antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunofluorescent analysis of HepG2 cells, using PDE12 antibody Catalog No:113656 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

Ig Type

Rabbit IgG


Applications

ELISA, WB, IF, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000



IgG (red).

IP Result of anti-PDE12 (IP:Catalog No:113656,

PDE12 antibody



Catalog Number: 113656

4ug; Detection:Catalog No:113656 1:500) with mouse brain tissue lysate 4000ug.