

PRODH Polyclonal Antibody



Catalog Number: 164656

Product name

PRODH Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PRODH

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 321-600 of human PRODH (NP_057419.4).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - PRODH Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PRODH antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.




Western blot - PRODH Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PRODH antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.

PRODH Polyclonal Antibody



Catalog Number: 164656



Immunohistochemistry - PRODH Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using PRODH antibody at dilution of 1:100 (40x lens).