

# RHEB antibody

Catalog Number: 114656

## Product name

RHEB antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RHEB Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, DU 145 cells, fetal human brain tissue, HeLa cells. Positive IHC detected in mouse brain tissue. Observed molecular weight by Western-blot: 20 kDa

## Preparation

This antibody was obtained by immunization of RHEB recombinant protein (Accession Number: NM\_005614). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114656(RHEB antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:114656(RHEB Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:114656(RHEB Antibody) at dilution of 1:50

# RHEB antibody



Catalog Number: 114656

---

(under 40x lens)