

# Anti-IL31 antibody



Catalog Number: 176656

## Product name

Anti-IL31 antibody

## Specificity

Human

## Antibody description

Mouse monoclonal antibody to IL31

## Preparation

This antigen of this antibody was recombinant protein

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

WB, ELISA

## Dilutions

WB: 1:2,000

ELISA: 1:5,000-1:10,000

## Validations




Fig1: Western blot analysis on bioactive recombinant human IL-31 protein using anti-IL-31 Mouse mAb (Cat. # 176656#).