

# CYB5R1 antibody



Catalog Number: 109656

## Product name

CYB5R1 antibody

WB: 1:1000-1:10000

IP: 1:500-1:5000

## Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

## Antibody description

CYB5R1 Rabbit Polyclonal antibody. Positive IF detected in PC-3 cells. Positive IHC detected in human breast cancer tissue. Positive IP detected in human placenta tissue. Positive WB detected in human placenta tissue, human brain tissue, PC-3 cells. Observed molecular weight by Western-blot: 34 kDa

## Validations




## Preparation

This antibody was obtained by immunization of CYB5R1 recombinant protein (Accession Number: NM\_016243). Purification method: Antigen affinity purified.

human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:109656(CYB5R1 antibody) at dilution of 1:500

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:109656(CYB5R1 antibody) at dilution of 1:50

## Applications


ELISA, WB, IHC, IF, IP

## Dilutions


Recommended Dilution:

# CYB5R1 antibody

Catalog Number: 109656



(under 10x lens)



Immunofluorescent analysis of PC-3 cells using Catalog No:109656(CYB5R1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

IP Result of anti-CYB5R1 (IP:Catalog No:109656, 4ug; Detection:Catalog No:109656 1:1000) with human placenta tissue lysate 1520ug.