

LSS antibody

Catalog Number: 112356

Product name

LSS antibody

IHC: 1:20-1:200

IF: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

LSS Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive IF detected in HepG2 cells. Positive WB detected in HeLa cells, HepG2 cells, human testis tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 70-75 kDa




Preparation

This antibody was obtained by immunization of LSS recombinant protein (Accession Number: NM_001001438). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112356(LSS antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human liver using Catalog No:112356(LSS antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

Immunohistochemical of paraffin-embedded human liver using Catalog No:112356(LSS

LSS antibody

Catalog Number: 112356

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (-20oc Ethanol) fixed HepG2 cells using Catalog No:112356(LSS Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)