

# NUDT12 antibody



Catalog Number: 113356

## Product name

NUDT12 antibody

## Clonality

Polyclonal

## Immunogen

[Human NUDT12 Recombinant protein \(His tag\)](#)

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

NUDT12 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse liver tissue. Observed molecular weight by Western-blot: 45-50 kDa


## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of NUDT12 recombinant protein (Accession Number: XM\_005272095). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113356(NUDT12 Antibody) at dilution of 1:600