

Anti-FPGT-TNNI3K antibody



Catalog Number: 176356

Product name

Anti-FPGT-TNNI3K antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to FPGT-TNNI3K

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tnni3k. aa range:301-350

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

IHC-P: 1:100-1:300

Validations




Fig1:; Western blot analysis of Mouse-kidney lysis using TNNI3K antibody. Antibody was diluted at 1:500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000