

ILVBL antibody

Catalog Number: 111756

Product name

ILVBL antibody

Specificity

Human; other species not tested.

Antibody description

ILVBL Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, MCF-7 cells. Positive IHC detected in human heart tissue. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of ILVBL recombinant protein (Accession Number: XM_005259717). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:111756(ILVBL Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:111756(ILVBL Antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:111756(ILVBL Antibody) at dilution of 1:100 (under 40x lens)