

# VIPAS39 antibody



Catalog Number: 116756

## Product name

VIPAS39 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

VIPAS39 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, A431 cells, HEK-293 cells, HT-1080 cells, human kidney tissue, L02 cells. Positive IF detected in Hela cells, HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 51 kDa

## Preparation

This antibody was obtained by immunization of VIPAS39 recombinant protein (Accession Number: NM\_001193317). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116756(C14orf133 antibody) at dilution of 1:300



Immunofluorescent analysis of Hela cells, using C14orf133 antibody Catalog No:116756 at 1:25 dilution and FITC- labeled FITC-Donkey anti-Rabbit IgG (green).