

Anti-DKK3 antibody



Catalog Number: 102066

Product name

Anti-DKK3 antibody

Immunogen

[Mouse DKK3 \(His Tag\) recombinant protein](#)

Specificity

Mouse Dkk-3 / Dickkopf-3 / DKK3

Antibody description

Rabbit polyclonal to DKK3

Preparation

Produced in rabbits immunized with purified, recombinant Mouse Dkk3 / Dickkopf-3 (rM Dkk3; NP_056629.1; Met 1-Ile 349). Dkk3 specific IgG was purified by mouse Dkk3 affinity chromatography

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.5-1.0 µg/mL

This antibody can be used at 0.5-1.0 µg/mL with the appropriate secondary reagents to detect Mouse Dkk-3. The detection limit for Mouse Dkk-3 is 0.00245 ng/well.

IHC-P: 0.1-1 µg/mL

Validations



DKK-3 / Dkk3 / Dickkopf-3 Antibody, Rabbit PAb, Antigen Affinity Purified

Immunochemical staining of mouse DKK3 in mouse brain with rabbit polyclonal antibody (0.2 µg/mL, formalin-fixed paraffin embedded sections).