

General Information

Gene Name Synonym

C-reactive protein(1-205)

Protein Construction

A DNA sequence encoding the human CRP (P02741-1) (Met1-Pro224) was expressed and purified.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 19

Molecular Mass

The recombinant human CRP consists of 206

amino acids and has a predicted molecular mass of 23 kDa. As a result of glycosylation, the apparent molecular mass of CRP is approximately 26 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 0.2M NaCl, 5mM CaCl₂, pH 8.0

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human CRP / C-Reactive Protein SDS-PAGE