

# ZMAT2 antibody



Catalog Number: 117066

## Product name

ZMAT2 antibody

## Immunogen

[Human ZMAT2 Recombinant protein \(His tag\)](#)

## Specificity

Human, Rat; other species not tested.

## Antibody description

ZMAT2 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, C6 cells, HeLa cells, Jurkat cells, MCF-7 cells. Observed molecular weight by Western-blot: 24 kDa

## Preparation

This antibody was obtained by immunization of ZMAT2 recombinant protein (Accession Number: NM\_144723). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:117066(ZMAT2 Antibody) at dilution of 1:600