

# Anti-p38 (5A1) Mouse antibody



Catalog Number: 168066

## Product name

Anti-p38 (5A1) Mouse antibody

## Specificity

Human, Mouse, Rat, Monkey

## Antibody description

p38 (5A1) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human P38 MAPK  
Pprotein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB

## Dilutions

WB: 1/500

## Validations



Western blot analysis of extracts from C2C12,PC12,HepG2,COS7,C6,3T3 and HeLa cell lysates using p38 MAPK mouse mAb (1:500 diluted).Predicted band size:40KDa.Observed band size:40KDa.