

Anti-p38 (5A1) Mouse antibody



Catalog Number: 168066

Product name

Anti-p38 (5A1) Mouse antibody

Ig Type

Mouse IgG1

Specificity

Human, Mouse, Rat, Monkey

Applications

WB

Antibody description

p38 (5A1) Mouse monoclonal antibody

Dilutions

WB: 1/500

Preparation

Antigen: Purified recombinant human P38 MAPK
Pprotein fragments expressed in E.coli.

Validations

Formulation


Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal



Western blot analysis of extracts from C2C12,PC12,HepG2,COS7,C6,3T3 and HeLa cell lysates using p38 MAPK mouse mAb (1:500 diluted).Predicted band size:40KDa.Observed band size:40KDa.