

Anti-PECAM1 antibody



Catalog Number: 142466

Product name

Anti-PECAM1 antibody

Immunogen

[Human CD31/PECAM-1 \(His Tag\) recombinant protein](#)

Specificity

Human CD31 / PECAM-1

Antibody description

Rabbit Monoclonal to Human PECAM1

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human CD31 / PECAM-1 .

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

ICC/IF, IF

Dilutions

ICC/IF: 10-25 µg/mL

Validations



CD31 / PECAM-1 Antibody, Rabbit MAb, Immunofluorescence