

Anti-LAIR2 antibody

Catalog Number: 176466

Product name

Anti-LAIR2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to LAIR2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human lair2. aa range:21-70

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western blot analysis of RAT-KIDNEY lysis using LAIR2 antibody. Antibody was diluted at 1:1000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000